




RGS14 Antibody

Product Code	CSB-PA019646ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43566
Immunogen	Recombinant Human Regulator of G-protein signaling 14 protein (230-510AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of G-protein signaling 14 (RGS14), RGS14
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGS14

Image



Western blot

All lanes: RGS14 antibody at 3.88μg/ml

Lane 1: 293T whole cell lysate

Lane 2: HeLa whole cell lysate

Lane 3: Jurkat whole cell lysate

Lane 4: SH-SY5Y whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 62, 22, 45, 37 kDa

Observed band size: 62 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.