





RHOBTB1 Antibody

Product Code	CSB-PA019684ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94844
Immunogen	Recombinant Human Rho-related BTB domain-containing protein 1 protein (240-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related BTB domain-containing protein 1, RHOBTB1, KIAA0740
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHOBTB1


Image



Western blot
 All lanes: RHOBTB1 antibody IgG at 3.17µg/ml +
 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 79 kDa
 Observed band size: 79 kDa



Immunohistochemistry of paraffin-embedded
 human thyroid tissue using CSB-
 PA019684ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA019684ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.