




ADCY4 Antibody

Product Code	CSB-PA822782ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFM4
Immunogen	Recombinant Human Adenylate cyclase type 4 protein (808-1077AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenylate cyclase type 4 (EC 4.6.1.1) (ATP pyrophosphate-lyase 4) (Adenylate cyclase type IV) (Adenyl cyclase 4), ADCY4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADCY4

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA822782ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.