





CIR1 Antibody

Product Code	CSB-PA773051DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86X95
Immunogen	Recombinant Human Corepressor interacting with RBPJ 1 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Corepressor interacting with RBPJ 1 (CBF1-interacting corepressor) (Recepin), CIR1, CIR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CIR1

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA773051DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA773051DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.