





MESP1 Antibody

Product Code	CSB-PA871559ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRJ9
Immunogen	Recombinant Human Mesoderm posterior protein 1 protein (149-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mesoderm posterior protein 1 (Class C basic helix-loop-helix protein 5) (bHLHc5), MESP1, BHLHC5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MESP1

Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA871559ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA871559ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.