




NUMB Antibody

Product Code	CSB-PA016186DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49757
Immunogen	Recombinant Human Protein numb homolog protein (532-651AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein numb homolog (h-Numb) (Protein S171), NUMB, C14orf41
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NUMB

Image



Western blot

All lanes: NUMB antibody at 0.93µg/ml

Lane 1: Mouse brain tissue

Lane 2: A431 whole cell lysate

Lane 3: HepG2 whole cell lysate

Lane 4: Hela whole cell lysate

Lane 5: Mouse kidney tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 71, 66, 70, 65, 56, 54, 51,

49, 38 kDa

Observed band size: 71 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.