





RGS9 Antibody

Product Code	CSB-PA019662ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75916
Immunogen	Recombinant Human Regulator of G-protein signaling 9 protein (230-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of G-protein signaling 9 (RGS9), RGS9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGS9

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA019662ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA019662ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.