




ZIC3 Antibody

Product Code	CSB-PA026486ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60481
Immunogen	Recombinant Human Zinc finger protein ZIC 3 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein ZIC 3 (Zinc finger protein 203) (Zinc finger protein of the cerebellum 3), ZIC3, ZNF203
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZIC3

Image



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA026486ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.