





RPIA Antibody

Product Code	CSB-PA020104ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49247
Immunogen	Recombinant Human Ribose-5-phosphate isomerase protein (82-311AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribose-5-phosphate isomerase (EC 5.3.1.6) (Phosphoriboisomerase), RPIA, RPI
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RPIA


Image



Western blot
 All lanes: RPIA antibody IgG at 3.29µg/ml +
 Mouse kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 33 kDa
 Observed band size: 33 kDa



Immunohistochemistry of paraffin-embedded
 human placenta tissue using CSB-
 PA020104ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA020104ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.