




# SNX18 Antibody

<b>Product Code</b>	CSB-PA850423ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96RF0
<b>Immunogen</b>	Recombinant Human Sorting nexin-18 protein (449-628AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Sorting nexin-18 (SH3 and PX domain-containing protein 3B) (Sorting nexin-associated Golgi protein 1), SNX18, SH3PXD3B SNAG1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SNX18

## Image



### Western blot

All lanes: SNX18 antibody at 1.6µg/ml

Lane 1: A549 whole cell lysate

Lane 2: Hela whole cell lysate


Lane 3: Mouse kidney tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 69, 70, 65 kDa


Observed band size: 69 kDa



### Immunohistochemistry of paraffin-embedded

human tonsil tissue using CSB-

PA850423ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA850423ESR1HU at dilution of 1:100