





CAMSAP2 Antibody

Product Code	CSB-PA804903ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08AD1
Immunogen	Recombinant Human Calmodulin-regulated spectrin-associated protein 2 protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calmodulin-regulated spectrin-associated protein 2 (Calmodulin-regulated spectrin-associated protein 1-like protein 1), CAMSAP2, CAMSAP1L1 KIAA1078
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAMSAP2

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA804903ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA804903ESR2HU at dilution of 1:100