





CORO1C Antibody

Product Code	CSB-PA891968DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULV4
Immunogen	Recombinant Human Coronin-1C protein (375-474AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coronin-1C (Coronin-3) (hCRNN4), CORO1C, CRN2 CRNN4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CORO1C

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA891968DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA891968DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.