





DNAJB2 Antibody

Product Code	CSB-PA007008DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25686
Immunogen	Recombinant Human DnaJ homolog subfamily B member 2 protein (215-324AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily B member 2 (Heat shock 40 kDa protein 3) (Heat shock protein J1) (HSJ-1), DNAJB2, HSJ1 HSPF3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJB2

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA007008DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA007008DSR1HU at dilution of 1:100