




EDIL3 Antibody

Product Code	CSB-PA007399ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43854
Immunogen	Recombinant Human EGF-like repeat and discoidin I-like domain-containing protein 3 protein (60-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	EGF-like repeat and discoidin I-like domain-containing protein 3 (Developmentally-regulated endothelial cell locus 1 protein) (Integrin-binding protein DEL1), EDIL3, DEL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EDIL3

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA007399ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.