





PIEZO1 Antibody

Product Code	CSB-PA852881ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92508
Immunogen	Recombinant Human Piezo-type mechanosensitive ion channel component 1 protein (2230-2420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Piezo-type mechanosensitive ion channel component 1 (Membrane protein induced by beta-amyloid treatment) (Mib) (Protein FAM38A), PIEZO1, FAM38A KIAA0233
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PIEZO1

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA852881ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA852881ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.