





FHL3 Antibody

Product Code	CSB-PA618796ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13643
Immunogen	Recombinant Human Four and a half LIM domains protein 3 protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Four and a half LIM domains protein 3 (FHL-3) (Skeletal muscle LIM-protein 2) (SLIM-2), FHL3, SLIM2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	FHL3


Image



Western blot
 All lanes: FHL3 antibody at 1.42µg/ml
 Lane 1: Human placenta tissue
 Lane 2: K562 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 31 kDa
 Observed band size: 31 kDa



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA618796ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA618796ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.