




SNRPB2 Antibody

Product Code	CSB-PA362460ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08579
Immunogen	Recombinant Human U2 small nuclear ribonucleoprotein B" protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	U2 small nuclear ribonucleoprotein B" (U2 snRNP B"), SNRPB2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SNRPB2

Image



Western blot
 All lanes: SNRPB2 antibody IgG at 2.26µg/ml +
 Mouse kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 25 kDa
 Observed band size: 25 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.