




SNRPB2 Antibody

Product Code	CSB-PA362460ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08579
Immunogen	Recombinant Human U2 small nuclear ribonucleoprotein B" protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	U2 small nuclear ribonucleoprotein B" (U2 snRNP B"), SNRPB2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SNRPB2

Image



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA362460ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.