





ERVFRD-1 Antibody

Product Code	CSB-PA352689ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60508
Immunogen	Recombinant Human Syncytin-2 protein (16-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syncytin-2 (Endogenous retrovirus group FRD member 1) (Envelope polyprotein) (HERV-FRD) (HERV-FRD_6p24.1 provirus ancestral Env polyprotein) [Cleaved into: Surface protein (SU); Transmembrane protein (TM)], ERVFRD-1, ERVFRDE1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ERVFRD-1

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA352689ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA352689ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.