





NRBP1 Antibody

Product Code	CSB-PA883410ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHY1
Immunogen	Recombinant Human Nuclear receptor-binding protein (406-535AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor-binding protein, NRBP1, BCON3 NRBP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NRBP1

Image



Western blot
 All lanes: NRBP1 antibody at 1.99µg/ml +
 HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 60 kDa
 Observed band size: 60 kDa



Immunohistochemistry of paraffin-embedded
 human testis tissue using CSB-
 PA883410ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.