





CHRNA6 Antibody

Product Code	CSB-PA620896ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15825
Immunogen	Recombinant Human Neuronal acetylcholine receptor subunit alpha-6 protein (26-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuronal acetylcholine receptor subunit alpha-6, CHRNA6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CHRNA6

Image



Western blot
 All lanes: CHRNA6 antibody at 2.29µg/ml
 Lane 1: A431 whole cell lysate
 Lane 2: THP-1 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 57, 56 kDa
 Observed band size: 57 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA620896ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.