




ITGAX Antibody

Product Code	CSB-PA011879ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20702
Immunogen	Recombinant Human Integrin alpha-X protein (100-350AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-X (CD11 antigen-like family member C) (Leu M5) (Leukocyte adhesion glycoprotein p150,95 alpha chain) (Leukocyte adhesion receptor p150,95) (CD antigen CD11c), ITGAX, CD11C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGAX

Image



Western blot

All lanes: ITGAX antibody IgG at 1.59µg/ml +


Mouse kidney tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 128 kDa


Observed band size: 128 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA011879ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA011879ESR1HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA011879ESR1HU at dilution of 1: 100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.