




THBD Antibody

Product Code	CSB-PA023486DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07204
Immunogen	Recombinant Human Thrombomodulin protein (200-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombomodulin (TM) (Fetomodulin) (CD antigen CD141), THBD, THRM
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	THBD

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA023486DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.