




CD276 Antibody

Product Code	CSB-PA745871ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5ZPR3
Immunogen	Recombinant Human CD276 antigen protein (29-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD276 antigen (4Ig-B7-H3) (B7 homolog 3) (B7-H3) (Costimulatory molecule) (CD antigen CD276), CD276, B7H3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD276

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA745871ESR2HU at dilution of 1:100