




# DSC2 Antibody

<b>Product Code</b>	CSB-PA007190ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q02487
<b>Immunogen</b>	Recombinant Human Desmocollin-2 protein (520-690AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Desmocollin-2 (Cadherin family member 2) (Desmocollin-3) (Desmosomal glycoprotein II) (Desmosomal glycoprotein III), DSC2, CDHF2 DSC3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	DSC2

## Image



Western blot  
 All lanes: DSC2 antibody at 3.46µg/ml  
 Lane 1: Rat heart tissue  
 Lane 2: A549 whole cell lysate  
 Lane 3: HepG2 whole cell lysate  
 Lane 4: Jurkat whole cell lysate  
 Lane 5: 293T whole cell lysate  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 100, 94 kDa  
 Observed band size: 94 kDa

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.