




DSC2 Antibody

Product Code	CSB-PA007190ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02487
Immunogen	Recombinant Human Desmocollin-2 protein (520-690AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Desmocollin-2 (Cadherin family member 2) (Desmocollin-3) (Desmosomal glycoprotein II) (Desmosomal glycoprotein III), DSC2, CDHF2 DSC3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DSC2

Image



Western blot

All lanes: DSC2 antibody at 3.46µg/ml

Lane 1: Rat heart tissue

Lane 2: A549 whole cell lysate

Lane 3: HepG2 whole cell lysate

Lane 4: Jurkat whole cell lysate

Lane 5: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 100, 94 kDa

Observed band size: 94 kDa