




DSC2 Antibody

Product Code	CSB-PA007190ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02487
Immunogen	Recombinant Human Desmocollin-2 protein (520-690AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Desmocollin-2 (Cadherin family member 2) (Desmocollin-3) (Desmosomal glycoprotein II) (Desmosomal glycoprotein III), DSC2, CDHF2 DSC3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DSC2

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA007190ESR2HU at dilution of 1:100