




# HHIP Antibody

<b>Product Code</b>	CSB-PA850419ESR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96QV1
<b>Immunogen</b>	Recombinant Human Hedgehog-interacting protein (18-220AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Rat
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Hedgehog-interacting protein (HHIP) (HIP), HHIP, HIP
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Stem Cells
<b>Target Names</b>	HHIP

## Image



### Western blot


All lanes: HHIPantibody IgG at 1.92μg/ml + Rat heart tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 79, 37 kDa

Observed band size: 79 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA850419ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA850419ESR2HU at dilution of 1:100