



HHIP Antibody


Product Code	CSB-PA850419ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QV1
Immunogen	Recombinant Human Hedgehog-interacting protein (18-220AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hedgehog-interacting protein (HHIP) (HIP), HHIP, HIP
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	HHIP

Image




Western blot

All lanes: HHIP antibody IgG at 1.92µg/ml + Rat heart tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 79, 37 kDa
 Observed band size: 79 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA850419ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA850419ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.