




SLC22A8 Antibody

Product Code	CSB-PA855047ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCC7
Immunogen	Recombinant Human Solute carrier family 22 member 8 protein (473-542AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 22 member 8 (Organic anion transporter 3) (hOAT3), SLC22A8, OAT3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC22A8

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA855047ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.