




TRHDE Antibody

Product Code	CSB-PA890740ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKU6
Immunogen	Recombinant fusion protein of human TRHDE
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:3000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 50% glycerol, pH7.4.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRHDE

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate(30µg), SY5Y whole cell lysate(30µg), Hela whole cell lysate(30µg), PC-3 whole cell lysate(30µg), THP-1 whole cell lysate(30µg), Mouse Heart tissue lysate(30µg)

All lanes: TRHDE antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 117, 123 kDa

Observed band size: 117 kDa

Exposure time: 120s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.