





ADAM19 Antibody

Product Code	CSB-PA880929ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H013
Immunogen	Recombinant Human Disintegrin and metalloproteinase domain-containing protein 19 protein (203-460AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disintegrin and metalloproteinase domain-containing protein 19 (ADAM 19) (EC 3.4.24.-) (Meltrin-beta) (Metalloprotease and disintegrin dendritic antigen marker) (MADDAM), ADAM19, MLTNB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADAM19

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA880929ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA880929ESR1HU at dilution of 1:100