





ASIC5 Antibody

Product Code	CSB-PA865174ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY37
Immunogen	Recombinant Human Acid-sensing ion channel 5 protein (140-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acid-sensing ion channel 5 (ASIC5) (Amiloride-sensitive cation channel 5) (Human intestine Na(+) channel) (HINaC), ASIC5, ACCN5 HINAC
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ASIC5

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA865174ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA865174ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.