





CD53 Antibody

Product Code	CSB-PA004944ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19397
Immunogen	Recombinant Human Leukocyte surface antigen CD53 protein (100-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte surface antigen CD53 (Cell surface glycoprotein CD53) (Tetraspanin-25) (Tspan-25) (CD antigen CD53), CD53, MOX44 TSPAN25
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD53

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA004944ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA004944ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.