





DSC1 Antibody

Product Code	CSB-PA007189ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08554
Immunogen	Recombinant Human Desmocollin-1 protein (480-690AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Desmocollin-1 (Cadherin family member 1) (Desmosomal glycoprotein 2/3) (DG2/DG3), DSC1, CDHF1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DSC1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA007189ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA007189ESR1HU at dilution of 1:100