




DUOX2 Antibody

Product Code	CSB-PA007229ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRD8
Immunogen	Recombinant Human Dual oxidase 2 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual oxidase 2 (EC 1.11.1.-) (EC 1.6.3.1) (Large NOX 2) (Long NOX 2) (NADH/NADPH thyroid oxidase p138-tox) (NADPH oxidase/peroxidase DUOX2) (NADPH thyroid oxidase 2) (Thyroid oxidase 2) (p138 thyroid oxidase), DUOX2, LNOX2 THOX2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DUOX2

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA007229ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.