




EFNB3 Antibody

Product Code	CSB-PA007467ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15768
Immunogen	Recombinant Human Ephrin-B3 protein (28-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin-B3 (EPH-related receptor transmembrane ligand ELK-L3) (EPH-related receptor tyrosine kinase ligand 8) (LERK-8), EFNB3, EPLG8 LERK8
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EFNB3

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA007467ESR2HU at dilution of 1:100