





FCGR2C Antibody

Product Code	CSB-PA008542ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31995
Immunogen	Recombinant Human Low affinity immunoglobulin gamma Fc region receptor II-c protein (43-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin gamma Fc region receptor II-c (IgG Fc receptor II-c) (CDw32) (Fc-gamma RII-c) (Fc-gamma-RIIc) (FcRII-c) (CD antigen CD32), FCGR2C, CD32 FCG2 IGFR2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCGR2C

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA008542ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA008542ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.