



GRIN3B Antibody

Product Code	CSB-PA009917ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60391
Immunogen	Recombinant Human Glutamate receptor ionotropic, NMDA 3B protein (180-410AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor ionotropic, NMDA 3B (GluN3B) (N-methyl-D-aspartate receptor subtype 3B) (NMDAR3B) (NR3B), GRIN3B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIN3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.