




# KIR3DL3 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA850818ESR2HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q8N743  |
| <b>Immunogen</b>           | Recombinant Human Killer cell immunoglobulin-like receptor 3DL3 protein (26-322AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Tested Applications</b> | ELISA, IHC; Recommended dilution: IHC:1:20-1:200  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.   |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Killer cell immunoglobulin-like receptor 3DL3 (CD158 antigen-like family member Z) (Killer cell inhibitory receptor 1) (CD antigen CD158z), KIR3DL3, CD158Z KIR3DL7 KIRC1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Immunology  |
| <b>Target Names</b>        | KIR3DL3   |

## Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA850818ESR2HU at dilution of 1:100