





ADGRL2 Antibody

Product Code	CSB-PA013060ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95490
Immunogen	Recombinant Human Adhesion G protein-coupled receptor L2 protein (400-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adhesion G protein-coupled receptor L2 (Calcium-independent alpha-latrotoxin receptor 2) (CIRL-2) (Letrophilin homolog 1) (Letrophilin-2) (Lectomedin-1), ADGRL2, KIAA0786 LEC1 LPHH1 LPHN2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ADGRL2

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA013060ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA013060ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.