




LRP4 Antibody

Product Code	CSB-PA013099ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75096
Immunogen	Recombinant Human Low-density lipoprotein receptor-related protein 4 protein (1500-1725AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low-density lipoprotein receptor-related protein 4 (LRP-4) (Multiple epidermal growth factor-like domains 7), LRP4, KIAA0816 LRP10 MEGF7
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LRP4

Image



Western blot


All lanes: LRP4 antibody at 1.39µg/ml + PC-3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 213 kDa

Observed band size: 213 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA013099ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.