




# NOX3 Antibody

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA015960ESR1HU                                                                                   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                        |
| <b>Uniprot No.</b>         | Q9HBY0                                                                                               |
| <b>Immunogen</b>           | Recombinant Human NADPH oxidase 3 protein (223-395AA)                                                |
| <b>Raised In</b>           | Rabbit                                                                                               |
| <b>Species Reactivity</b>  | Human                                                                                                |
| <b>Tested Applications</b> | ELISA, IHC; Recommended dilution: IHC:1:20-1:200                                                     |
| <b>Form</b>                | Liquid                                                                                               |
| <b>Conjugate</b>           | Non-conjugated                                                                                       |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                                    |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                            |
| <b>Isotype</b>             | IgG                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                           |
| <b>Alias</b>               | NADPH oxidase 3 (EC 1.6.3.-) (Mitogenic oxidase 2) (MOX-2) (gp91phox homolog 3) (GP91-3), NOX3, MOX2 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                 |
| <b>Research Area</b>       | Neuroscience                                                                                         |
| <b>Target Names</b>        | NOX3                                                                                                 |

## Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA015960ESR1HU at dilution of 1:100

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.