




NRXN1 Antibody

Product Code	CSB-PA347832ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULB1
Immunogen	Recombinant Human Neurexin-1 protein (31-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurexin-1 (Neurexin I-alpha) (Neurexin-1-alpha), NRXN1, KIAA0578
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NRXN1

Image



Western blot

All lanes: NRXN1 antibody at 1.57µg/ml

Lane 1: Jurkat whole cell lysate

Lane 2: Hela whole cell lysate

Lane 3: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 162, 165, 170, 47, 16, 51 kDa

Observed band size: 162 kDa