



SCN1A Antibody

Product Code	CSB-PA020834ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35498
Immunogen	Recombinant Human Sodium channel protein type 1 subunit alpha protein (274-399AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium channel protein type 1 subunit alpha (Sodium channel protein brain I subunit alpha) (Sodium channel protein type I subunit alpha) (Voltage-gated sodium channel subunit alpha Nav1.1), SCN1A, NAC1 SCN1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SCN1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.