





# SLCO2B1 Antibody

<b>Product Code</b>	CSB-PA021746ESR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O94956
<b>Immunogen</b>	Recombinant Human Solute carrier organic anion transporter family member 2B1 protein (130-210AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Solute carrier organic anion transporter family member 2B1 (Organic anion transporter B) (OATP-B) (Organic anion transporter polypeptide-related protein 2) (OATP-RP2) (OATPRP2) (Solute carrier family 21 member 9), SLCO2B1, KIAA0880 OATP2B1 OATPB SLC21A9
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SLCO2B1

## Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA021746ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA021746ESR2HU at dilution of 1:100

---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.