





SLCO2B1 Antibody

Product Code	CSB-PA021746ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94956
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 2B1 protein (130-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 2B1 (Organic anion transporter B) (OATP-B) (Organic anion transporter polypeptide-related protein 2) (OATP-RP2) (OATPRP2) (Solute carrier family 21 member 9), SLCO2B1, KIAA0880 OATP2B1 OATPB SLC21A9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLCO2B1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA021746ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA021746ESR2HU at dilution of 1:100