




TAS1R3 Antibody

Product Code	CSB-PA752534ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7RTX0
Immunogen	Recombinant Human Taste receptor type 1 member 3 protein (400-570AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 1 member 3 (Sweet taste receptor T1R3), TAS1R3, T1R3 TR3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS1R3

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA752534ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.