




TAS1R3 Antibody

Product Code	CSB-PA752534ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7RTX0
Immunogen	Recombinant Human Taste receptor type 1 member 3 protein (400-570AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS1R3

Image



Western Blot

Positive WB detected in: PC-3 whole cell lysate(20µg), MCF7 whole cell lysate(20µg), JK whole cell lysate(20µg)

All lanes: TAS1R3 antibody at 1:500


Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 94 kDa

Observed band size: 100 kDa

Exposure time: 300s



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA752534ESR2HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.