




TAS1R3 Antibody

Product Code	CSB-PA752534ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7RTX0
Immunogen	Recombinant Human Taste receptor type 1 member 3 protein (400-570AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 1 member 3 (Sweet taste receptor T1R3), TAS1R3, T1R3 TR3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS1R3

Image



Western blot

All lanes: TAS1R3 antibody at 1.48µg/ml

Lane 1: K562 whole cell lysate

Lane 2: HL60 whole cell lysate

Lane 3: Hela whole cell lysate


Lane 4: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 94 kDa

Observed band size: 94 kDa



Immunohistochemistry of paraffin-embedded

human pancreatic cancer using CSB-

PA752534ESR2HU at dilution of 1:100

